ICP - Mass Spectrometry

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Ruth Merrified¹Jamie Lead¹Chady Stephan²

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1. Sigg, L.; Behra, R.; Groh, K.; Isaacson, C.; Odzak, N.; Piccapietra, F.; Rohder, L.; Schug, H.; Yue, Y.; Schirmer, K., Chemical Aspects of Nanoparticle Ecotoxicology. *Chimia* 2014, 68, (11), 806-811.
2. Zhao, L. L.; Kim, T. H.; Kim, H. W.; Ahn, J. C.; Kim, S. Y., Enhanced cellular uptake and phototoxicity of Verteporfin-conjugated gold nanoparticles as theranostic nanocarriers for targeted photodynamic therapy and imaging of cancers. *Materials Science & Engineering C-Materials for Biological Applications* 2016, 67, 611-622.
3. Akrami, M.; Balalaie, S.; Hosseinkhani, S.; Alipour, M.; Salehi, F.; Bahador, A.; Haririan, I., Tuning the anticancer activity of a novel pro-apoptotic peptide using gold nanoparticle platforms. *Scientific Reports* 2016, 6, 12.
4. Egger, A. E.; Rappel, C.; Jakupec, M. A.; Hartinger, C. G.; Heffeter, P.; Keppler, B. K., Development of an experimental protocol for uptake studies of metal compounds in adherent tumor cells. *Journal of Analytical Atomic Spectrometry* 2009, 24, (1), 51-61.
5. Journal of Physics: Conference Series 522, (2014), 012058.
6. Taylor, N. S.; Merrifield, R.; Williams, T. D.; Chipman, J. K.; Lead, J. R.; Viant, M. R., Molecular toxicity of cerium oxide nanoparticles to the freshwater alga *Chlamydomonas reinhardtii* is associated with supra-environmental exposure concentrations. *Nanotoxicology* 2016, 10, (1), 32-41.
7. Wolfbeis, O. S., An overview of nanoparticles commonly used in fluorescent bioimaging. *Chemical Society Reviews* 2015, 44, (14), 4743-4768.
8. Hineman, A.; Stephan, C., Effect of dwell time on single particle inductively coupled plasma mass spectrometry data acquisition quality. *Journal of Analytical Atomic Spectrometry* 2014, 29, (7), 1252-1257.
9. Montano, M. D.; Olesik, J. W.; Barber, A. G.; Challis, K.; Ranville, J. F., Single Particle ICP-MS: Advances toward routine analysis of nanomaterials. *Analytical and Bioanalytical Chemistry* 2016, 408, (19), 5053-5074.
10. "A Comparison of Microsecond vs. Millisecond Dwell Times on Particle Number Concentration Measurements by Single Particle ICP-MS", PerkinElmer application note, 2016.
11. "Syngistix Nano Application Module for Single Particle ICP-MS", Product Note, PerkinElmer Inc., 2014.
12. Abad-Álvo, I., Peña Vazquez, E., Bolea, E., Bermejo-Barrera, P., Castillo, J., Laborda, F. *Anal. Bioanal. Chem.*, 2016, 408, 5089.

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©2017, PerkinElmer, Inc. 版权所有。保留所有权利。PerkinElmer® 是 PerkinElmer, Inc. 的注册商标。所有其他商标均为其各自所有者的财产。所有解释权归PerkinElmer。