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


<sup>1</sup>Department of Chemistry and Environmental  
 Science Centre, Brunel University of Science  
 and Technology

<sup>2</sup>Centre for Single Nanoparticles, Single Cell, and  
 Single Molecules Microscopy (CSM3)

<sup>3</sup>Department of Civil and Environmental  
 Engineering, Utah State University

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Effectiveness of Three Water Treatment Plants Removing TiO<sub>2</sub> Particles and Dissolved Ti.

Plant	Pre/ Post Treatment	Most Frequent Size (nm)	Particle Concentration (particles/mL)	Dissolved Concentration (µg/L)
1	Pre Post	170 < MDL	432,000 < MDL	17.9 1.21
2	Pre Post	156 < MDL	451,000 < MDL	11.7 1.17
3	Pre Post	153 76	425,000 17,237	10.6 < MDL

Sample	Au			Ag			TiO <sub>2</sub>		
	Most Freq Size (nm)	Part Conc. Spike Recovery	Diss Conc. Spike Recovery	Most Freq Size (nm)	Part Conc. Spike Recovery	Diss Conc. Spike Recovery	Most Freq Size (nm)	Part Conc. Spike Recovery	Diss Conc. Spike Recovery
1	98	97%	80%	98	97%	80%	102	9%	84%
2	97	88%	84%	97	88%	84%	87	6%	88%
3	101	94%	89%	101	94%	89%	87	6%	112%

# Tracking Nanoparticles Dissolution/Aggregation

APPLICATION NOTE

ICP - Mass Spectrometry

Authors:  
 Denise Mitrano  
 James F. Kazelle  
 Department of Chemistry and Geochemistry  
 Colorado School of Mines  
 Golden, CO USA

Chady Stepien  
 PerkinElmer, Inc.  
 Shelton, CT

### Quantitative Evaluation of Nanoparticle Dissolution Kinetics using Single Particle ICP-MS: A Case Study with Silver Nanoparticles

**Introduction**  
 Accurate data on engineered nanoparticle (ENP) environmental behavior and the interplay between ENP size, dissolution rate, agglomeration, and interaction with the sample matrix is critical to appropriately characterize the risks these novel materials may pose to environmental health. The advancement of the single particle ICP-MS (SP-ICP-MS) technique is a great benefit for the study of ENPs in natural systems at environmentally relevant (ng/L) concentrations. Previous studies may have obscured environmentally-relevant transformations because of artificially high ENP concentrations used in the experiments<sup>1</sup>. Therefore, the SP-ICP-MS method is at the forefront to garner the type of information most relevant for environmental risk assessments, namely the precise tracking of changes in ENP size, associated dissolved metal concentration, and determining polydispersity of an ENP sample, all at dilute concentrations in complex solutions. Because dissolution rate is surface-area controlled, the time to complete dissolution is highly dependent on the initial and (potentially stable) intermediate particle sizes. By measuring the change in particle size, as well as the evolution of Ag(Iaq) in solution, using SP-ICP-MS, potential pitfalls related to loss of Ag<sup>+</sup> to experimental materials and to other environmental surfaces, such as suspended sediments or biota in the case of complex matrices, may be avoided.


Figure 2. Particle size distribution of Ag ENP suspended in various water chemistries (DI, tap, and creek water) over 24 h. Evidence of decreasing particle diameter with time through particle oxidation and dissolution in some samples (e.g. DI and tap water) with less change in particle size observed in other samples (e.g. creek water).

# Silver Nanoparticles in Wastewaters

## APPLICATION NOTE

### ICP – Mass Spectrometry

#### Authors:

Mehnoosh Anadi

Subhas Ghoshal

McGill University

Montreal, Canada

Chady Stephan

PerkinElmer, Inc.

Shelton, CT

## Measurement and Analysis of Silver Nanoparticles in Wastewaters with Single Particle ICP-MS

### Introduction

The drastic increase in production and consumption of engineered nanoparticles (ENPs) has raised the concern and questions about their release into the environment and

potential harm to aquatic and terrestrial species. The characteristic properties of nanoparticles, such as small size and high specific surface area and reactivity, make them desirable for their use in various products.

Silver (Ag) nanoparticles are among the most commonly used nanoparticles in consumer products due to their antimicrobial properties. Therefore, it is expected that Ag ENPs will find their way into the environment, necessitating a way to accurately and rapidly detect and characterize them in a variety of environmental matrices. Work has already been performed demonstrating the ability to successfully characterize Ag ENPs in a variety of water samples<sup>1</sup> and biological media which may be exposed to Ag ENPs in the environment<sup>2</sup>.



Figure 2. Measured Ag particle size distribution in effluent wastewater diluted 1000 times.

Table 2. Results from Analysis of Wastewater Samples Spiked with Ag NPs.

Sample	Mean Particle Size (nm)	Spiked Particle Conc. (particles/mL)	Measured Particle Conc. (particles/mL)
Effluent Wastewater	66.3±0.2	50,000	54,691±1185
Mixed Liquor Wastewater	63.7±0.4	50,000	53,123±1216

# Bioavailability of Nanoparticles in Tomato Plants

## Gold Nanoparticle Uptake by Tomato Plants Characterized by Single Particle ICP-MS

### Introduction

With the increasing use of engineered nanoparticles (ENPs) in a variety of products and processes, there is concern about the release of ENPs into and impact on the environment.

One aspect of the environmental concern is the potential for ENPs to make their way to plants via migration through water and/or soil. If ENPs end up in food crops, this is a potential pathway to human exposure.

The challenge arises in how to measure ENPs in plant materials and, more specifically, in sample preparation. To our knowledge, current sample preparation techniques have limited capability to conserve the concentration and characteristics of nanoparticles (NPs) once they enter plant tissues, as they mainly depend on acid digestion. These limitations can be avoided by careful choice of the ENP extraction procedure and performing the analysis with single particle ICP-MS (SP-ICP-MS), the combination of which will preserve the particle size information, allow for rapid analysis of a large number of samples, and yield results on the particle size, concentration, and size distribution.


Table 1. NexION 300/350D Instrumental and Analytical Parameters.


Parameter	Value
Nebulizer	Concentric (glass)
Nebulizer Flow	1.08 L/min
Spray Chamber	Baffled Cyclonic (glass)
ICP RF Power	1600 W
Analyte	Au
Mass	197 amu
Dwell Time	0.1 ms
Settling Time	0 ms
Sampling Time	100 sec
Number of Data Points Acquired	1 million per sample
Au Density	19.3 g/cm <sup>3</sup>

# Metal containing Nanoparticles in Sludges

Environmental Risk Implications of Metals in Sludges from Waste Water Treatment Plants: The Discovery of Vast Stores of Metal-Containing Nanoparticles

Ferun Tom,<sup>1</sup> Yi Yang,<sup>1,2\*</sup> Jiguan Feng,<sup>1</sup> Zhenhui Niu,<sup>1</sup> Hui Pan,<sup>1</sup> Yukun Qin,<sup>1</sup> Xingyun Gao,<sup>1</sup> Xiangzhou Meng,<sup>3</sup> Min Liu,<sup>1</sup> and Michael F. Hochella<sup>1,4</sup>


# Nanoparticles in Sunscreens


Sunscreen	SPF	TiO <sub>2</sub> Content (%)
1	60+	4.9
2	50	6
3	45	0
4	50	6
5	45+	5.1

Sunscreen	Dilution Factor	Most Frequent Size (nm)	Particle Size Distribution (nm)	Particle Concentration (particles/mL)
1	100,000	32	24 – 58	102,229
2	100,000	34	24 – 64	117,252
3*	100,000	-	-	-
4	100,000	33	24 – 61	63,000
5	20,000	42	28 – 67	198,462


Sample	Analyte	Most Freq. Size (nm)	Mean Size (nm)	No. of Peaks	Mean Inten. (counts)	Part. Conc. (parts/mL)
Ti48 for 2_rep1 20000x	Ti 47.948	37	67	11239	104.07	266258
Ti48 for 2_rep2 20000x	Ti 47.948	37	65	132943	88.76	431565
Ti48 for 2_rep3 20000x	Ti 47.948	37	63	12658	80.83	412501
Ti48 for 4_rep1 20000x	Ti 47.948	37	67	7279	116.26	237209
Ti48 for 4_rep2 20000x	Ti 47.948	37	65	8609	108.3	280551
Ti48 for 4_rep3 20000x	Ti 47.948	36	65	8442	100.53	275109
Ti48 for 5_rep1 20000x	Ti 47.948	39	72	4092	154.89	133351
Ti48 for 5_rep2 20000x	Ti 47.948	38	67	6895	95.01	234695
Ti48 for 5_rep3 20000x	Ti 47.948	38	68	6000	100.8	198462
QC 5 ppb	Ti 47.948	75	136	9	493.22	293
Ti48 for 1_rep1 100000x	Ti 47.948	25	43	2969	32.86	96784
Ti48 for 1_rep2 100000x	Ti 47.948	27	47	3137	33.48	102229
Ti48 for 1_rep3 100000x	Ti 47.948	27	48	3653	37.94	119045
Ti48 for 2_rep1 100000x	Ti 47.948	37	64	2758	83.73	89678
Ti48 for 2_rep2 100000x	Ti 47.948	33	58	3611	34.90	117076
Ti48 for 2_rep3 100000x	Ti 47.948	32	58	3598	88.1	117252
Ti48 for 4_rep1 100000x	Ti 47.948	32	59	1923	77.61	62667
Ti48 for 4_rep2 100000x	Ti 47.948	30	59	2577	74.62	83580
Ti48 for 4_rep3 100000x	Ti 47.948	37	66	1863	104.89	60712
Ti48 for 5_rep1 10000x	Ti 47.948	39	70	10376	103.89	338135
Ti48 for 5_rep2 10000x	Ti 47.948	39	68	13160	95.08	428860
Ti48 for 5_rep3 10000x	Ti 47.948	37	64	17141	80.72	958393
QC 5 ppb	Ti 47.948	75	109	12	306.58	391

# Nanoparticles in Biological Samples



<p>Urine</p> 	<p>1/10 dilution in 0.5% <math>\text{NH}_4\text{OH}</math> + 0.1% Triton</p>
<p>Blood</p> 	<p>1/20 dilution in 0.5% <math>\text{NH}_4\text{OH}</math> + 0.1% Triton</p>
<p>Tissues</p> 	<p>24-hour digestion with Tetramethylammonium Hydroxide (TMAH -20% w/w) followed by centrifugation</p>
<p>Organisms</p> 	<p>24-hour digestion with Tetramethylammonium Hydroxide (TMAH - 20% w/w) followed by centrifugation</p>

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
## New Research Evaluating Cisplatin Uptake in Ovarian Cancer Cells by Single Cell ICP-MS

Recent cancer chemotherapy drugs in the ovarian cancer category effectiveness is due to its ability to bind to the DNA, resulting in DNA-platinum (Pt) adducts, which bend the DNA. The cells must then repair the DNA damage, otherwise DNA replication is inhibited, resulting in cell death.

Many cancers are initially sensitive to platinum-based treatment, but patients frequently relapse with tumors displaying resistance to further platinum therapy. Cisplatin drug resistance is due to three major molecular mechanisms: increased DNA repair, increased cellular inactivation, and altered cellular accumulation. Increased cellular uptake or increased cellular export of cisplatin contributes the mechanisms involved in altered cellular accumulation.

### Introduction

Cisplatin, carboplatin, and oxaliplatin are the most widely used of platinum-based cancer chemotherapy drugs in the ovarian cancer category. Cisplatin effectiveness is due to its ability to bind to the DNA, resulting in DNA-platinum (Pt) adducts, which bend the DNA. The cells must then repair the DNA damage, otherwise DNA replication is inhibited, resulting in cell death.

Many cancers are initially sensitive to platinum-based treatment, but patients frequently relapse with tumors displaying resistance to further platinum therapy. Cisplatin drug resistance is due to three major molecular mechanisms: increased DNA repair, increased cellular inactivation, and altered cellular accumulation. Increased cellular uptake or increased cellular export of cisplatin contributes the mechanisms involved in altered cellular accumulation.

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## Quantification of Gold Nanoparticle Uptake into Cancer Cells using Single Cell ICP-MS

### Introduction

Cancer is a disease that is characterized by the uncontrolled growth and spread of abnormal cells. According to the American Cancer Society, cancer is the second most common cause of death in the US.<sup>1</sup> Current treatments for various cancers include surgery, radiation, immunotherapy, and chemotherapy. Although these conventional therapies may improve patients' overall survival and quality of life, they also have several limitations. For example, in conventional cancer chemotherapy, small-molecule-based cancer therapeutics distribute non-specifically throughout the entire human body. The consequence is that these drugs do not only kill cancer cells but also destroy healthy cells in the body causing severe side effects for cancer patients,<sup>2</sup> leading to a need for new therapies that can target diseased cells.

Nanoparticles (NPs) have been of significant interest over the last two decades as they offer attractive benefits for drug delivery to overcome limitations in conventional chemotherapy.<sup>3</sup> Nanoparticles can be engineered to carry both drugs and imaging probes to simultaneously detect and treat cancer. They may also be designed to specifically target diseased tissues and cells in the body. A number of nanoparticle-based cancer therapeutics have been approved for clinical use and/or are currently under development.<sup>4,5</sup> Advantages that rationally engineered nanoparticles may offer over conventional small-molecule drugs include: (i) prolonged circulation time in the body; (ii) reduction of nonspecific cellular uptake and undesirable off-target and side effects; and (iii) improvement in cellular interactions through specific cancer cell targeting moieties.


Figure 5. SC-ICP-MS analysis of human bladder cancer cells (T24) incubated with PEGylated 50-nm (gold core size) AuNPs for 4 h, using a PerkinElmer NexION system. (A) Histogram of gold mass distribution for individual T24 cells – blue curves highlight populations of cells containing one or more AuNPs per cell. (B) Normalized frequency distribution histogram of cells containing one or more AuNPs (norm. to number of cells containing one AuNP; n=3) – bars indicate average values with standard deviations.

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